



Horizontal Field Mapping System

FMU-MRI Product Introduction



RESONANCE



RESEARCH

MAGNETIC FIELD MAPPING SYSTEM FOR HORIZONTAL BORE MAGNETS

- RRI's FMU-MRI system, used for mapping static magnetic fields in Horizontal-bore, MRI magnets
- Operates on ^1H NMR and a heli-spherical NMR probe actuator.
- Digitally synthesized RF with FAST-SWEEP NMR Spectrometer with 10 MHz range
- The spectrometer operates from below 5 MHz (0.12 Tesla) to above 1500 MHz (35 Tesla).

MAGNETIC FIELD MAPPING SYSTEM

MAIN OPERATING PROPERTIES

- Operation at proton (H-1) resonance frequency
- Equipped with a heli-spherical path electro-mechanical probe and designed for a single sample operation.
- High resolution stepper motor under computer control with 1600 steps per revolution
- The heli-cylindrical probe actuator was designed for use in bore diameters between 200 and 1200 mm.
- Use of the spherical-path probe is very efficient and result in data suitable for axial and radial component analysis simultaneously.
- System completes a full mapping and analysis cycle in less than 30 minutes.

Mapping Probe Parameters

Older version shown



Probe trajectory
HELISPHERICAL


Radius **75-250mm**

Pitch **9.8 deg**

Length **1600 mm**

Sample size **1 mm³**

HELISPHERICAL SAMPLE TRAJECTORY



Mapper Principles of Operation

- Magnet center and probe orientation by difference plot, typically ZX
- Shim calibration by difference plot and spherical harmonic analysis
- Results stored in .STR files
- Shim convergence by matrix inversion using the .STR files
- Option for remote shim power supply operation
- Residual gradient calculation with SD error

Shim Calibration by Difference plot


- Acquire magnet Map 1 with shim setting to Value 1
- Acquire magnet Map 2 with shim setting to Value 2
- Automated difference map plot
- Spherical harmonic coefficient calculation
- Store data

Mapping setup of 9.4T/800

Courtesy Dr. Keith Thulborn UIC



Difference Plot of ZX Shim



Typical Spherical Harmonic Coefficient Analysis of a Map (shown 750 MHz vertical bore)

Results of Analysis

2007-12-24 14:59:04

c:\user\mit75sc\map16.dat

Gradient	z0 =	-1478.71	+/-	3.55	Hz/cm^n
Gradient	z =	-147.46	+/-	3.11	Hz/cm^n
Gradient	z2 =	151.04	+/-	4.17	Hz/cm^n
Gradient	z3 =	-21.79	+/-	4.91	Hz/cm^n
Gradient	z4 =	-3.02	+/-	6.64	Hz/cm^n
Gradient	z5 =	0.0000			
Gradient	z6 =	0.0000			
Gradient	z7 =	0.0000			
Gradient	z8 =	0.0000			
Gradient	z9 =	0.0000			
Gradient	z10 =	0.0000			
Gradient	z11 =	0.0000			
Gradient	x =	-45.99	+/-	5.36	Hz/cm^n
Gradient	y =	93.22	+/-	5.38	Hz/cm^n
Gradient	zx =	878.85	+/-	8.01	Hz/cm^n
Gradient	zy =	-113.18	+/-	8.11	Hz/cm^n
Gradient	c2 =	4.64	+/-	6.61	Hz/cm^n
Gradient	s2 =	9.18	+/-	6.60	Hz/cm^n
Gradient	z2x =	3.55	+/-	6.08	Hz/cm^n
Gradient	z2y =	-5.07	+/-	6.29	Hz/cm^n
Gradient	zc2 =	-5.68	+/-	5.58	Hz/cm^n
Gradient	zs2 =	-3.81	+/-	5.56	Hz/cm^n
Gradient	c3 =	-1.90	+/-	5.66	Hz/cm^n
Gradient	s3 =	5.52	+/-	5.64	Hz/cm^n
Gradient	z3x =	-6.89	+/-	8.04	Hz/cm^n
Gradient	z3y =	4.97	+/-	8.58	Hz/cm^n
Gradient	z2c2 =	-27.59	+/-	7.30	Hz/cm^n
Gradient	z2s2 =	-11.60	+/-	7.27	Hz/cm^n
Gradient	zc3 =	2.31	+/-	6.58	Hz/cm^n
Gradient	zs3 =	5.50	+/-	6.50	Hz/cm^n

Sum of squared residuals is 44306.7

Correlation coefficient is 0.999517


RMS uncertainty is 6.307




PRESS <Enter> TO CONTINUE

7 T 900 mm Bore Magnet – raw field


Courtesy of Dr. Alan Koretsky NIH




Plot of Final SC Convergence on 45 cm DSV




Final RT Convergence with 17 channel shims on 45 cm DSV




Final RT Convergence ppm contour plot (MathLab post-processing)




Final RT Convergence 1H Hz contour plot in X-Y plane



Final RT Convergence 1H Hz contour plot in Z-X plane



Final RT Convergence 1H Hz contour plot in Z-Y plane



Other Functions

- **Drift measurement**
 - Single point frequency plotted against time
- **Multiple map repeat and analysis**
- **Probe tuning**
- **Extended analysis and map editing options**
- **Shim simulations**
- **Shim power calculations**
- **Passive shimming as service by RRI**

Magnet Ramp Power Supply



MXR-800

- Output: 880A @ 12V
- Programmable ramp sequences
- Ramp rates adjustable up to 120A/min in 1A increments
- 20mA current resolution
- Quench detection circuit
- Programmable switch heater current up to 1A
- Reinforced shippable enclosure

Thank you for your attention

